TESTING PLATFORM	PRINCIPLE	ADVANTAGES	DISADVANTAGES
Lateral flow immunoassays: SARS-CoV-2 antigen assays	<ol> <li>The individual's sample is placed in a reagent tube, which contains a lysis buffer that disrupts the cells to expose the nucleocapsid protein inside.</li> <li>When this solution containing patient specimen and reagent is added to the testing cassette, the sample moves along a test strip.</li> <li>A colour change is occurs when a positive sample reaches the part of the strip where the antibody against the N protein is located.</li> </ol>	<ul> <li>Low cost.</li> <li>Fast, the results are obtained in less than 30 minutes.</li> <li>Easy to perform.</li> <li>Doesn't require complex equipment.</li> </ul>	<ul> <li>The sensitivity and specificity depends on the brand used.</li> <li>They are used as screening tests and PCR may be needed to confirm the results.</li> </ul>
Lateral flow immunoassays: SARS-CoV-2 antibody assays	<ol> <li>The test principle is based on the receptor-binding domain (RBD) of the spike and nucleocapsid proteins (antigens).</li> <li>The cassette has dye pad sections containing Recombinant 2019-novel coronavirus nucleocapsid protein and or Recombinant 2019-novel coronavirus Spike Protein (Si Subunit).</li> <li>The blood/serum/plasma sample is placed in the testing cassette, the sample moves along a test strip.</li> <li>A colour change occurs when a positive sample (a positive sample is one containing antibodies against SARS-COV-2 virus) reaches the part of the strip where the SARS-CoV-2 antigen is located.</li> </ol>	<ul> <li>Low cost.</li> <li>Fast, the results are obtained in less than 30 minutes</li> <li>Doesn't require complex equipment</li> </ul>	<ul> <li>Determines exposure but does not predict the current infection.</li> <li>The sensitivity and specificity depends on the brand used.</li> </ul>
Reverse transcriptase Real- time PCR assays	<ol> <li>RNA is extracted from the sample usually using a kit based method.</li> <li>Reverse transcriptase Real-time PCR is then performed on the extracted nucleic acids.</li> <li>SARS-CoV-2 gene targets commonly used:</li> <li>nucleocapsid, envelope protein, RNA dependent RNA polymerase, the viral open reading frames.</li> <li>The Taqman, FRET and LAMP chemistries are usually used.</li> <li>Automated and manual PCR methods are available.</li> </ol>	<ul> <li>Gold standard for SARS-CoV- 2 testing therefore it serves as a confirmatory test for other tests.</li> <li>High specificity and sensitivity.</li> </ul>	<ul> <li>Requires expensive material and equipment.</li> <li>Needs expertise and training</li> <li>Long turn around (1-4hours ).</li> <li>Higher cost.</li> <li>Requires a functional laboratory.</li> </ul>